

# Make your routine work more comfortable

## SAVE 15% ON YOUR FAVOURITE EPPENDORF CENTRIFUGE

Extraordinary versatility centrifuge for centrifugation of tubes and bottles from 0.2 ml to 750 ml, but also all types of MTP, PCR, cell culture, or Deepwell Plates.

- For centrifugation of up to 3L sample
- High centrifugation speed of up to  $20913 \times g$  (14000 rpm) depending on the rotor equipped
- Ergonomic user-friendly and easy: one key dedicated to each standard function
- Display of the speed in rpm or in  $\times g$
- User can adjust the centrifuging radius, for more precision
- "At set rpm" function starts the timer when the selected rotational speed is reached for reproducible centrifugation runs
- Acceleration and braking ramps independently adjustable (10 steps each)
- Emergency stop switch on front side of the unit
- Automatic rotor recognition and imbalance detection for maximum operational safety
- Soft-touch lid closure, allows closing the centrifuge with ease
- 3 year warranty



### 5810 centrifuge:

- Noise level: <65 dB(A)
- Dimensions (W × D × H): 54 × 61 × 35 cm
- Weight: 68 kg

### 5810R centrifuge:

- Refrigerated rotor, temperature settings from  $-9\text{ }^{\circ}\text{C}$  to  $40\text{ }^{\circ}\text{C}$
- FastTemp function for fast pre-cooling, maintains constant  $4\text{ }^{\circ}\text{C}$ , even at maximum speed
- Continuous cooling maintains a constant temperature after the run has been completed - your samples stay cool
- Noise level: <59 dB(A)
- Dimensions (W × D × H): 70 × 61 × 35 cm
- Weight: 99 kg

**Please use campaign code when ordering: CM0430**

**Turn to see more information ►►►**

# Make your routine work more comfortable

Item no.	Description	EUR list price	EUR discounted price
5810000320	Centrifuge 5810 G, 230 V/50-60 Hz, incl. rotor A-4-62 and 15/50ml adapters	7.505	6.379
5810000420	Centrifuge 5810 G, 230 V/50-60 Hz, incl. rotor A-4-81 and 15/50ml adapters	8.437	7.172
5810000620	Centrifuge 5810 G, 230 V/50-60 Hz, with rotor S-4-104 incl. adapter for round-bottom tubes diameter 13 / 16 mm	9.036	7.681
5811000325	Centrifuge 5810 R G (refrigerated), 230 V/50-60 Hz, incl. rotor A-4-62 and 15/50ml adapters	10.743	9.132
5811000420	Centrifuge 5810 R G (refrigerated), 230 V/50-60 Hz, incl. rotor A-4-81 and 15/50ml adapters	11.660	9.911
5811000465	Centrifuge 5810 R G (refrigerated), 230 V/50-60 Hz, incl. rotor A-4-81 and 15/50ml adapters, with GB-plug	11.660	9.911

**Please use campaign code when ordering: CM0430**

## Rotor A-4-62:

- Speed up to 3220 x g (4000 rpm)
- Accommodates tubes and flasks from 0.2 ml to 250 ml as well as plates



## Rotor A-4-81:

- Speed up to 3220 x g (4000 rpm)
- Accommodates tubes and flasks from 0.2 ml to 500 ml as well as plates
- Mixed loading with 2 plate buckets and 2 rectangular tube buckets



## Rotor S-4-104:

- Speed up to 3214 x g (3900 rpm)
- Accommodates tubes and bottles from 0.2 ml to 750 ml as well as plates
- Round bucket design for up to 36 % energy savings
- Mixed loading with 2 plate buckets and 2 round tube bucket



Note: Aerosol-tight caps optionally available for the buckets, please see Cat. No. 5820747005, 5820748001 and 5810710006.